Mini Mold™
TECHNIQUE

Mold-A-Lateral™

5 Steps - 5 Minutes

Turns a peg/small lateral
into a full-size lateral

Lateral Bonding Kit

8 Tips total, 2 Handles and
Detailed Step-By-Step Instructions per kit:

4 Right Lateral Bonder Tips

- 2 Large (8.0mm width)
- 2 Small (7.0mm width)

4 Left Lateral Bonder Tips

- 2 Large (8.0mm width)
- 2 Small (7.0mm width)

REF 51388

ORTHO TECHNOLOGY®

17401 Commerce Park Blvd. • Tampa, Florida, 33647 USA

(813) 991-5896 • Fax: (813) 991-5986

Secure Online Ordering: www.orthotechnology.com

Mini-Mold™ Technique

5 Minute Lateral Bonding Kit

Reusable Molds

Instructions for Use

Preliminary Set-up

If the tooth is in contact, open at least 0.5mm interproximal space. If a normal size lateral exists on the opposite side, measure the width and select the appropriate small (7.0mm) or large (8.0mm) mold. If bilateral bonding is necessary, measure the central incisor and use the mold that is approximately 2.0mm narrower. If the clinical crown is exceedingly short, trimming of the gingival portion of the mold may be necessary.

Bonding Technique

Insert the Lateral Bonder Tip into the handle. Estimate the amount of resin to place in the mold by assessing the size and volume of the tooth to be bonded. Example: for a peg tooth, fill about ½ of the mold. Fill less for a larger tooth.

Helpful Tip

For extra-large laterals, a slight squeeze of the large mold with the thumb and forefinger before curing will slightly widen the width of the mold.

Step 1 - Prep

Etch and prime (seal) the tooth.

Step 2 - Test Fit

Select the tip size, fill with composite and test fit the Mold-A-Lateral until incisal edge aligns with centrals. Withdraw the mold, remove any gingival excess if the mold is overfilled or add more composite if mold is underfilled.

Step 3 - Cure

Place the mold back onto the tooth in the desired position and cure with curing light (10 to 20 seconds depending on the light intensity and material).

Step 4 - Remove

When the Light-Cure Resin is fully cured, slowly remove the tip of the mold by pulling the mold straight off of the tooth.

Step 5 - Finish

Use a finishing bur for any contour adjustments. Interproximal stripping may be utilized to narrow the width of the restoration, if necessary.

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