**Intensiv SA**  
Via al Molino 107  
6926 Montagnola  
Switzerland  
www.intensiv.ch

## Intensiv Ortho-Strips System

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Intensiv Ortho-Strips System, Double-Sided, Ref. OS08OP-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS80XC-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS60C-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS40M-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS25F-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS15POL-DS

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Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-R, rechts
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-R, rechts
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-R, rechts
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-R, rechts
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-R, rechts
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-L, links
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-R, rechts

**Set**

Intensiv Ortho-Strips Double-Sided Set, Ref. 060B (OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS)
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15POL-L, OS15POL-R, OS25F-L, OS25F-R, OS40M-L, OS40M-R)
Intensiv Ortho-Strips One-Sided Set, Ref. 065 (OS80XC-R, OS80XC-L, OS60C-R, OS60C-L, OS40M-R, OS40M-L, OS15POL-R, OS15POL-L, OS08OP-R, OS08OP-L)
Intensiv Ortho-Strips Set IPR, Ref. 070 (OS80XC-DS, OS60C-DS, OS40M-DS, OS15POL-DS, OS08OP-DS)

Neue Ref.	Vorherige Ref.
OS80XC-DS	OS80C
OS80XC-R	OS80C-R
OS80XC-L	OS80C-L
OS08OP-DS	OS08OP
OS08OP-R	OS08OP-R
OS08OP-L	OS08OP-L
OS15POL-DS	OS15
OS15POL-R	OS15R
OS15POL-L	OS15L
OS25F-DS	OS25
OS25F-R	OS25R
OS25F-L	OS25L
OS40M-DS	OS40
OS40M-R	OS40R
OS40M-L	OS40L
OS60C-DS	OS60
OS60C-R	OS60R
OS60C-L	OS60L

## Anwendungshinweise Instructions for use Mode d’emploi

## Intensiv Ortho-Strips System

University of Zurich

Oscillating instruments for easy and rapid opening of tight interdental con- tact points and for the unilateral and bilateral interproximal reduction of the tooth size and finishing and polishing of interproximal surfaces in Or- thodontics.

To be used in combination with the contra-angles:

- Intensiv Swingle, WG-69 LT Ortho, with light
- Intensiv Swingle, WG-69 A Ortho, without light

Alternatively:

- KaVo head 61 LG
- KaVo contra-angle E61C

**Product description**

- Double-sided diamond-coated oscillating metal strip:
  - Intensiv Ortho-Strips Opener OS08OP-DS: serrated and diamond coated on both sides, 8µm
  - Intensiv Ortho-Strips Extra Coarse OS80XC-DS: perforated and diamond coated on both sides, 80µm
  - Intensiv Ortho-Strips Coarse OS60C-DS: diamond coated on both sides, 60µm
  - Intensiv Ortho-Strips Medium OS40M-DS: diamond coated on both sides, 40µm
  - Intensiv Ortho-Strips Fine OS25F-DS: diamond coated on both sides, 25µm
  - Intensiv Ortho-Strips Polisher OS15POL-DS: diamond coated on both sides, 15µm
- One-sided diamond-coated oscillating metal strip (L= coated on left side. R= coated on right side):
  - Intensiv Ortho-Strips Opener, OS08OP-L or OS08OP-R: serrated and di- amond coated on one side, 8µm
  - Intensiv Ortho-Strips Extra Coarse, OS80XC-L or OS80XC-R: perforated and diamond coated on one side, 80µm
  - Intensiv Ortho-Strips Coarse, OS60C-L or OS60C-R: diamond coated on one side, 60µm
  - Intensiv Ortho-Strips Medium, OS40M-L or OS40M-R: diamond coated on one side, 40µm
  - Intensiv Ortho-Strips Fine, OS25F-L or OS25F-R: diamond coated on one side, 25µm
  - Intensiv Ortho-Strips Polisher, OS15POL-L or OS15POLR: diamond coated on one side, 15µm
- Strip height: 3.7 mm
- Length of the diamond-coated part: 13 mm
- Sterilizable and reusable

**Indications**

- Unilateral and bilateral initial opening of tight interdental contact points at the beginning of the IPR orthodontics protocol
- Unilateral and bilateral interproximal enamel reduction during the IPR or- thodontics protocol
- Unilateral and bilateral finishing and polishing of the interproximal sur- faces.

**Recommendations for use**

- Apply Intensiv Ortho-Strips Opener at the beginning of the IPR orthodon- tics protocol for the opening of tight interproximal contact points.
- Subsequently apply Intensiv Ortho-Strips according to the IPR protocol for orthodontical treatments.
- In case of very tight contact points separate the teeth by placing a wooden wedge cervically. This makes stripping more comfortable both for operator and patient.
- Set the micro motor at 40'000 RPM/min max. With Intensiv Swingle this corresponds to 20'000 oscillations per minute approx. (2:1).
- Activate water spray during the entire treatment (minimum 50 ml/min).
- Regulate water flow in case of water excess or shortage.
- Insert the strip in a cervical direction while the contra-angle is moving. This prevents the instruments from jamming.
- Apply a load not greater than 100 grams, up to a maximum bending of the strip of approximately 1 mm.
- Do not stop the strip in the interdental space.
- Remove the strip from the interdental space whilst the contraangle is mov- ing.

**Maintenance and Sterilization**

- Instruments are delivered non sterile. They must be sterilized prior to first use and disinfected, cleaned, and sterilized after each use.
- Clean the instruments and remove debris after each use, so as to maintain their abrasive properties.
- Use a single-pulsed or fractionated vacuum autoclave and subvacuum dry- ing. Chemiclave sterilizers may also be used. Hot air sterilizers are not suit- able for diamond instruments.
- Sterilization must be carried out according to validated procedures.
- Concerning the sterilization process we refer to the ISO standard 17664. We suggest to follow the indications below:
  - Cycles at 134°C:
  - Tmin = 134°C – Tmax = 138°C
  - Pressure = 3.15 bar abs
  - Time = 4 min (raisable)
  - Cycles at 121°C:
  - Tmin = 121°C – Tmax = 125°C
  - Pressure = 2.10 bar abs
  - Time = 16 min (raisable)

**Risk Warnings**

- Avoid bending and levering: use of excessive force (strong levering effect with the contra-angle) could lead to the breakage of the instrument.
- Avoid jamming or levering actions when oscillating, as this increases the risk of instrument breakage.
- Having passed over the tooth equator line with Intensiv Ortho-Strips Open- er, pay attention not to injure the soft tissues.
- Oscillating over 20'000 oscillations per minute could cause breakage of the instrument.
- Avoid temperatures above 180°C which may affect the durability of the instrument.
- Avoid applying loads greater than 100 grams as this could cause damage to the instrument and/or the tooth.
- To ensure traceability of the instruments during their entire application, we recommend keeping the packaging (lot number).

**Packaging of 1, 3, 6 or 12 pieces**

Intensiv Ortho-Strips System, Double-Sided, Ref. OS08OP-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS80XC-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS60C-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS40M-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS25F-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS15POL-DS

**Packaging of 1, 3, 6 or 12 pieces**

Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-R, right
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-R, right
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-R, right
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-R, right
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-R, right
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-L, left
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-R, right

**Set**

Intensiv Ortho-Strips Double-Sided Set, Ref. 060B (OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS)
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15POL-L, OS15POL-R, OS25F-L, OS25F-R, OS40M-L, OS40M-R)
Intensiv Ortho-Strips One-Sided Set, Ref. 065 (OS80XC-R, OS80XC-L, OS60C-R, OS60C-L, OS40M-R, OS40M-L, OS15POL-R, OS15POL-L, OS08OP-R, OS08OP-L)
Intensiv Ortho-Strips Set IPR, Ref. 070 (OS80XC-DS, OS60C-DS, OS40M-DS, OS15POL-DS, OS08OP-DS)

New Ref.	Previous Ref.
OS80XC-DS	OS80C
OS80XC-R	OS80C-R
OS80XC-L	OS80C-L
OS08OP-DS	OS08OP
OS08OP-R	OS08OP-R
OS08OP-L	OS08OP-L
OS15POL-DS	OS15
OS15POL-R	OS15R
OS15POL-L	OS15L
OS25F-DS	OS25
OS25F-R	OS25R
OS25F-L	OS25L
OS40M-DS	OS40
OS40M-R	OS40R
OS40M-L	OS40L
OS60C-DS	OS60
OS60C-R	OS60R
OS60C-L	OS60L

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Intensiv Ortho-Strips System, Double-Sided, Ref. OS08OP-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS80XC-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS60C-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS40M-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS25F-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS15POL-DS

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Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-R, droite
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-R, droite
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-R, droite
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-R, droite
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-R, droite
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-L, gauche
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-R, droite

**Set**

Intensiv Ortho-Strips Double-Sided Set, Ref. 060B (OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS)
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15POL-L, OS15POL-R, OS25F-L, OS25F-R, OS40M-L, OS40M-R)
Intensiv Ortho-Strips One-Sided Set, Ref. 065 (OS80XC-R, OS80XC-L, OS60C-R, OS60C-L, OS40M-R, OS40M-L, OS15POL-R, OS15POL-L, OS08OP-R, OS08OP-L)
Intensiv Ortho-Strips Set IPR, Ref. 070 (OS80XC-DS, OS60C-DS, OS40M-DS, OS15POL-DS, OS08OP-DS)

Réf. Nouvelle	Réf. Précédente
OS80XC-DS	OS80C
OS80XC-R	OS80C-R
OS80XC-L	OS80C-L
OS08OP-DS	OS08OP
OS08OP-R	OS08OP-R
OS08OP-L	OS08OP-L
OS15POL-DS	OS15
OS15POL-R	OS15R
OS15POL-L	OS15L
OS25F-DS	OS25
OS25F-R	OS25R
OS25F-L	OS25L
OS40M-DS	OS40
OS40M-R	OS40R
OS40M-L	OS40L
OS60C-DS	OS60
OS60C-R	OS60R
OS60C-L	OS60L

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Intensiv Ortho-Strips System, Double-Sided, Ref. OS08OP-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS80XC-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS60C-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS40M-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS25F-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS15POL-DS

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Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-R, destra
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-R, destra
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-R, destra
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-R, destra
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-R, destra
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-L, sinistra
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-R, destra

#### Set

Intensiv Ortho-Strips Double-Sided Set, Ref. 060B (OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS)
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15POL-L, OS15POL-R, OS25F-L, OS25F-R, OS40M-L, OS40M-R)
Intensiv Ortho-Strips One-Sided Set, Ref. 065 (OS80XC-R, OS80XC-L, OS60C-R, OS60C-L, OS40M-R, OS40M-L, OS15POL-R, OS15POL-L, OS08OP-R, OS08OP-L)
Intensiv Ortho-Strips Set IPR, Ref. 070 (OS80XC-DS, OS60C-DS, OS40M-DS, OS15POL-DS, OS08OP-DS)

Ref. Nuova	Ref. Precedente
OS80XC-DS	OS80C
OS80XC-R	OS80C-R
OS80XC-L	OS80C-L
OS08OP-DS	OS08OP
OS08OP-R	OS08OP-R
OS08OP-L	OS08OP-L
OS15POL-DS	OS15
OS15POL-R	OS15R
OS15POL-L	OS15L
OS25F-DS	OS25
OS25F-R	OS25R
OS25F-L	OS25L
OS40M-DS	OS40
OS40M-R	OS40R
OS40M-L	OS40L
OS60C-DS	OS60
OS60C-R	OS60R
OS60C-L	OS60L

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#### Embalajes per 1, 3, 6 o 12 piezas

Intensiv Ortho-Strips System, Double-Sided, Ref. OS08OP-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS80XC-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS60C-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS40M-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS25F-DS
Intensiv Ortho-Strips System, Double-Sided, Ref. OS15POL-DS

#### Embalajes per 1, 3, 6 o 12 piezas

Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS08OP-R, derecha
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS80XC-R, derecha
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS60C-R, derecha
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS40M-R, derecha
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS25F-R, derecha
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-L, izquierda
Intensiv Ortho-Strips System, One-Sided, Ref. OS15POL-R, derecha

#### Set

Intensiv Ortho-Strips Double-Sided Set, Ref. 060B (OS60C-DS, OS40M-DS, OS25F-DS, OS15POL-DS)
Intensiv Ortho-Strips One-Sided Set, Ref. 064 (OS15POL-L, OS15POL-R, OS25F-L, OS25F-R, OS40M-L, OS40M-R)
Intensiv Ortho-Strips One-Sided Set, Ref. 065 (OS80XC-R, OS80XC-L, OS60C-R, OS60C-L, OS40M-R, OS40M-L, OS-15POL-R, OS15POL-L, OS08OP-R, OS08OP-L)
Intensiv Ortho-Strips Set IPR, Ref. 070 (OS80XC-DS, OS60C-DS, OS40M-DS, OS15POL-DS, OS08OP-DS)

Ref. Nueva	Ref. Anterior
OS80XC-DS	OS80C
OS80XC-R	OS80C-R
OS80XC-L	OS80C-L
OS08OP-DS	OS08OP
OS08OP-R	OS08OP-R
OS08OP-L	OS08OP-L
OS15POL-DS	OS15
OS15POL-R	OS15R
OS15POL-L	OS15L
OS25F-DS	OS25
OS25F-R	OS25R
OS25F-L	OS25L
OS40M-DS	OS40
OS40M-R	OS40R
OS40M-L	OS40L
OS60C-DS	OS60
OS60C-R	OS60R
OS60C-L	OS60L

 **Intensiv SA**  
Via al Molino 107  
6926 Montagnola  
Switzerland

Tel. +41 91 986 50 50  
Fax +41 91 986 50 59  
info@intensiv.ch  
www.intensiv.ch