

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Date of issue: 2020/05/07 Version: 1

#### **SECTION 1: Identification**

1.1. Identification

Product form : Mixture

Product name : Glass Lok Liquid

Product code : GLL15, GLL25, GLL15B, GLL25B

1.2. Recommended use and restrictions on use

Use of the substance/mixture : For Rx Only

#### 1.3. Supplier

#### Manufacturer:

Reliance Orthodontic Products, Inc. 1540 West Thorndale Ave.

Itasca, IL 60143 USA

630-773-4009, during normal business hours

www.RelianceOrthodontics.com

#### 1.4. Emergency telephone number

Emergency number : CHEMTREC - 24-Hour Hazmat Emergency Communications Center

Domestic: 1-800-424-9300 Outside the U.S.: 1-703-527-3887, collect calls accepted

# SECTION 2: Hazard(s) identification

#### 2.1. Classification of the substance or mixture

#### **GHS-US** classification

Skin corrosion/irritation, Category 1 H314 Causes severe skin burns and eye damage. Skin sensitisation, Category 1 H317 May cause an allergic skin reaction.

Full text of H statements : see section 16

#### 2.2. GHS Label elements, including precautionary statements

#### **GHS US labelling**

Hazard pictograms (GHS US)





Signal word (GHS US) : Danger

Hazard statements (GHS US) : H314 - Causes severe skin burns and eye damage.

H317 - May cause an allergic skin reaction.

Precautionary statements (GHS US) : P260 - Do not breathe vapours, gas.

P261 - Avoid breathing vapours.

P264 - Wash hands thoroughly after handling

P272 - Contaminated work clothing must not be allowed out of the workplace P280 - Wear protective gloves/protective clothing/eye protection/face protection

P301+P330+P331 - If swallowed: rinse mouth. Do NOT induce vomiting

P302+P352 - If on skin: Wash with plenty of water and soap

P303+P361+P353 - If on skin (or hair): Take off immediately all contaminated clothing. Rinse

skin with water/shower

P304+P340 - If inhaled: Remove person to fresh air and keep comfortable for breathing P305+P351+P338 - IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing. P310 - Immediately call Call a doctor, a POISON CENTER P321 - Specific treatment (see First aid measures on this label)

P333+P313 - If skin irritation or rash occurs: Get medical advice/attention.

P363 - Wash contaminated clothing before reuse.

P405 - Store locked up.

P501 - Dispose in a safe manner in accordance with local/national regulations

# 2.3. Other hazards which do not result in classification

No additional information available

### 2.4. Unknown acute toxicity (GHS US)

Not applicable

EN (English) Page 1

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

#### **SECTION 3: Composition/information on ingredients**

#### **Substances**

Not applicable

#### 3.2. **Mixtures**

Name	Product identifier	%	GHS-US classification
2-Hydroxyethyl Methacrylate	(CAS-No.) 868-77-9	10 - 30	Skin Irrit. 2, H315 Eye Irrit. 2, H319 Skin Sens. 1, H317
Glycerol DMA	(CAS-No.) 1830-78-0	10 - 30	Skin Irrit. 2, H315 Eye Irrit. 2A, H319 STOT SE 3, H335
10-Methacryloyloxydecyl Dihydrogen Phosphate	(CAS-No.) 85590-00-7	5 - 10	Skin Irrit. 2, H315 Eye Irrit. 2, H319 STOT SE 3, H335

Full text of hazard classes and H-statements : see section 16

#### **SECTION 4: First-aid measures**

#### **Description of first aid measures**

First-aid measures general Get medical advice/attention if you feel unwell. Never give anything by mouth to an

unconscious person. If medical advice is needed, have product container or label at hand.

Remove person to fresh air and keep comfortable for breathing. Get medical advice/attention if First-aid measures after inhalation

you feel unwell.

Wash immediately with lots of water. Soap may be used. If skin irritation or rash occurs: Get First-aid measures after skin contact medical advice/attention. Take off immediately all contaminated clothing and wash it before

If eye irritation persists: Get medical advice/attention. Rinse immediately with plenty of water. First-aid measures after eye contact

Obtain medical attention if pain, blinking or redness persists.

First-aid measures after ingestion If you feel unwell, seek medical advice. Rinse mouth with water. Do NOT induce vomiting.

#### Most important symptoms and effects (acute and delayed)

Potential adverse human health effects and

symptoms

: Based on available data, the classification criteria are not met.

: Causes eye irritation. May cause respiratory irritation. May cause an allergic skin reaction. Symptoms/effects

Symptoms/effects after inhalation : May cause respiratory irritation. May cause an allergic skin reaction.

Symptoms/effects after skin contact Causes skin irritation. Symptoms/effects after eye contact May cause eye irritation.

Symptoms/effects after ingestion : May cause irritation to the digestive tract.

#### Immediate medical attention and special treatment, if necessary

Treat symptomatically.

### **SECTION 5: Fire-fighting measures**

#### Suitable (and unsuitable) extinguishing media

Suitable extinguishing media : Dry chemical powder. Alcohol-resistant foam. Carbon dioxide. Water spray. Sand.

Unsuitable extinguishing media : Do not use a heavy water stream.

#### Specific hazards arising from the chemical

No additional information available

#### Special protective equipment and precautions for fire-fighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any

chemical fire. Prevent fire fighting water from entering the environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

#### **SECTION 6: Accidental release measures**

#### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Emergency procedures : Ventilate spillage area. Avoid contact with skin and eyes. Evacuate unnecessary personnel.

> EN (English) 2/7

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

#### 6.1.2. For emergency responders

Protective equipment : Do not attempt to take action without suitable protective equipment. For further information

refer to section 8: "Exposure controls/personal protection". Equip cleanup crew with proper

protection.

Emergency procedures : Ventilate area.

#### 6.2. Environmental precautions

Avoid release to the environment. Prevent entry to sewers and public waters. Notify authorities if liquid enters sewers or public waters.

#### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Take up liquid spill into absorbent material. Soak up spills with inert solids, such as clay or

diatomaceous earth as soon as possible. Collect spillage. Store away from other materials.

Other information : Dispose of materials or solid residues at an authorized site.

#### 6.4. Reference to other sections

For further information refer to section 13. See Heading 8. Exposure controls and personal protection.

#### SECTION 7: Handling and storage

#### 7.1. Precautions for safe handling

Precautions for safe handling : Wear personal protective equipment. Avoid contact with skin, eyes and clothing. Wash hands

and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Provide good ventilation in process area to prevent formation of vapour.

Avoid breathing fume/vapours.

Hygiene measures : Always wash hands after handling the product. Contaminated work clothing should not be

allowed out of the workplace. Wash contaminated clothing before reuse.

#### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep cool. Protect from sunlight. Keep container closed when not in use.

Incompatible products : Strong bases. Strong acids.
Incompatible materials : Sources of ignition. Direct sunlight.

#### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

# Glass Lok Liquid

No additional information available

# 2-Hydroxyethyl Methacrylate (868-77-9)

No additional information available

#### Glycerol DMA (1830-78-0)

No additional information available

#### 10-Methacryloyloxydecyl Dihydrogen Phosphate (85590-00-7)

No additional information available

#### 8.2. Appropriate engineering controls

Appropriate engineering controls : Ensure good ventilation of the work station.

Environmental exposure controls : Avoid release to the environment.

#### 8.3. Individual protection measures/Personal protective equipment

#### Personal protective equipment:

Avoid all unnecessary exposure.

#### Hand protection:

Protective gloves. Wear protective gloves/protective clothing/eye protection/face protection protective gloves

#### Eye protection:

Safety glasses. Chemical goggles or safety glasses

#### Skin and body protection:

Wear suitable protective clothing

EN (English) 3/7

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

#### Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Wear appropriate mask

#### Other information:

Do not eat, drink or smoke during use.

# **SECTION 9: Physical and chemical properties**

#### 9.1. Information on basic physical and chemical properties

Physical state: LiquidAppearance: Clear liquid.Colour: Blue, Light yellow

Odour : Acrylic

Odour threshold : No data available

pH : ≈ 1

Melting point : Not applicable
Freezing point : No data available
Boiling point : No data available
Flash point : No data available
Relative evaporation rate (butylacetate=1) : No data available

Flammability (solid, gas) : Not applicable. Non flammable.

Vapour pressure : No data available : No data available Relative vapour density at 20 °C : No data available Relative density Solubility : No data available No data available Log Pow Auto-ignition temperature : No data available Decomposition temperature : No data available Viscosity, kinematic : No data available : No data available Viscosity, dynamic **Explosive limits** No data available : No data available Explosive properties Oxidising properties : No data available

#### 9.2. Other information

No additional information available

### **SECTION 10: Stability and reactivity**

# 10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

### 10.2. Chemical stability

Stable under normal conditions of use. Not established.

### 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use. Not established.

### 10.4. Conditions to avoid

None under recommended storage and handling conditions (see section 7). Direct sunlight. Extremely high or low temperatures.

#### 10.5. Incompatible materials

Strong acids. Strong bases.

#### 10.6. Hazardous decomposition products

Under normal conditions of storage and use, hazardous decomposition products should not be produced. fume. Carbon monoxide. Carbon dioxide.

#### **SECTION 11: Toxicological information**

### 11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified

EN (English) 4/7

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Acute toxicity (dermal) : Not classified
Acute toxicity (inhalation) : Not classified

2-Hydroxyethyl Methacrylate (868-77-9)	
LD50 oral rat	5564 mg/kg bodyweight (Rat, Experimental value, Oral)
LD50 dermal rabbit > 5000 mg/kg (24 h, Rabbit, Male, Experimental value, Dermal)	
ATE US (oral)	5564 mg/kg bodyweight

# 10-Methacryloyloxydecyl Dihydrogen Phosphate (85590-00-7)

LD50 oral rat > mg/kg

Skin corrosion/irritation : Causes severe skin burns and eye damage.

pH: ≈ 1

Serious eye damage/irritation : Serious eye damage, category 1, implicit

pH: ≈ 1

Respiratory or skin sensitisation : May cause an allergic skin reaction.

Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : Not classified

### Glycerol DMA (1830-78-0)

STOT-single exposure May cause respiratory irritation.

# 10-Methacryloyloxydecyl Dihydrogen Phosphate (85590-00-7)

STOT-single exposure May cause respiratory irritation.

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified
Viscosity, kinematic : No data available

Potential adverse human health effects and

symptoms

: Based on available data, the classification criteria are not met.

Symptoms/effects : Causes eye irritation. May cause respiratory irritation. May cause an allergic skin reaction.

Symptoms/effects after inhalation : May cause respiratory irritation. May cause an allergic skin reaction.

Symptoms/effects after skin contact : Causes skin irritation.

Symptoms/effects after eye contact : May cause eye irritation.

Symptoms/effects after ingestion : May cause irritation to the digestive tract.

# **SECTION 12: Ecological information**

#### 12.1. Toxicity

Ecology - general : Before neutralisation, the product may represent a danger to aquatic organisms.

2-Hydroxyethyl Methacrylate (868-77-9)		
LC50 fish 1	> 100 mg/l (OECD 203: Fish, Acute Toxicity Test, 96 h, Oryzias latipes, Semi-static system, Fresh water, Experimental value, GLP)	
EC50 Daphnia 1	380 mg/l (OECD 202: Daphnia sp. Acute Immobilisation Test, 48 h, Daphnia magna, Static system, Fresh water, Experimental value, GLP)	
EC50 Daphnia 2	380 mg/l (EC50; OECD 202: Daphnia sp. Acute Immobilisation Test; 48 h; Daphnia magna; Static system; Fresh water; Experimental value)	
ErC50 (algae)	836 mg/l (OECD 201: Alga, Growth Inhibition Test, 72 h, Pseudokirchneriella subcapitata, Static system, Fresh water, Experimental value, GLP)	
Threshold limit algae 1	836 mg/l (ErC50; OECD 201: Alga, Growth Inhibition Test; 72 h; Pseudokirchneriella subcapitata; Static system; Fresh water; Experimental value)	
Threshold limit algae 2	345 mg/l (EbC50; OECD 201: Alga, Growth Inhibition Test; 72 h; Pseudokirchneriella subcapitata; Static system; Fresh water; Experimental value)	

EN (English) 5/7

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

10-Methacryloyloxydecyl Dihydrogen Phosphate (85590-00-7)	
EC50 Daphnia 1	> mg/l
NOEC chronic fish	48h 10 mg/l

# 12.2. Persistence and degradability

Blass Lok Liquid		
Persistence and degradability	Not established.	
2-Hydroxyethyl Methacrylate (868-77-9)		
Persistence and degradability Biodegradability in soil: no data available. Readily biodegradable in water.		
Glycerol DMA (1830-78-0)		
Persistence and degradability	Persistence and degradability Biodegradability in water: no data available.	

#### 12.3. Bioaccumulative potential

Glass Lok Liquid	
Bioaccumulative potential	Not established.
2-Hydroxyethyl Methacrylate (868-77-9)	
BCF fish 1	1.3 - 1.5 (Pisces, Calculated value)
Log Pow	0.42 (Experimental value, OECD 107: Partition Coefficient (n-octanol/water): Shake Flask Method, 25 °C)
Bioaccumulative potential	Not bioaccumulative.
Glycerol DMA (1830-78-0)	
Bioaccumulative potential	No bioaccumulation data available.

# 12.4. Mobility in soil

2-Hydroxyethyl Methacrylate (868-77-9)	
Ecology - soil	Adsorbs into the soil.

# 12.5. Other adverse effects

Other information : Avoid release to the environment.

# **SECTION 13: Disposal considerations**

#### 13.1. Disposal methods

Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions.

Product/Packaging disposal recommendations : Avoid release to the environment. Dispose in a safe manner in accordance with local/national

regulations.

Ecology - waste materials : Avoid release to the environment.

# SECTION 14: Transport information

#### **Department of Transportation (DOT)**

In accordance with DOT

Not regulated

### **Transportation of Dangerous Goods**

Not applicable

#### Transport by sea

Not applicable

#### Air transport

Not applicable

EN (English) 6/7

# Safety Data Sheet

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

#### SECTION 15: Regulatory information

#### 15.1. US Federal regulations

All components of this product are listed, or excluded from listing, on the United States Environmental Protection Agency Toxic Substances Control Act (TSCA) inventory except for:

١	10-Methacryloyloxydecyl Dihydrogen Phosphate	CAS-No. 85590-00-7	5 - 10%

This product or mixture is not known to contain a toxic chemical or chemicals in excess of the applicable de minimis concentration as specified in 40 CFR §372.38(a) subject to the reporting requirements of section 313 of Title III of the Superfund Amendments and Reauthorization Act of 1986 and 40 CFR Part 372.

#### 15.2. International regulations

#### **CANADA**

### 2-Hydroxyethyl Methacrylate (868-77-9)

Listed on the Canadian DSL (Domestic Substances List)

#### Glycerol DMA (1830-78-0)

Listed on the Canadian NDSL (Non-Domestic Substances List)

#### 10-Methacryloyloxydecyl Dihydrogen Phosphate (85590-00-7)

Not listed on the Canadian DSL (Domestic Substances List)/NDSL (Non-Domestic Substances List)

#### **EU-Regulations**

#### **National regulations**

No additional information available

#### 15.3. US State regulations

California Proposition 65 - This product does not contain any substances known to the state of California to cause cancer, developmental and/or reproductive harm

# **SECTION 16: Other information**

according to Federal Register / Vol. 77, No. 58 / Monday, March 26, 2012 / Rules and Regulations

Revision date : 09/30/2019 Other information : None.

#### Full text of H-statements:

H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.

#### SDS US (GHS HazCom 2012)

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product.

EN (English) 7/7